

Viscous Fluid Delivery

Viscous Fluid Cannula

High-pressure rated 4 mm beveled tip infusion cannula, ideal for injection of silicone oil and other viscous fluids, as well as routine saline infusion.

Flanged cannula designed to easily suture in place.

3 Cannulas per box, Sterile, Single Use

Reorder Number TR9849-51

Viscous Fluid Pack

Single-use procedure pack for injection and aspiration of silicone oil and other viscous fluids. Connects to Sonomed Escalon Viscous Fluid System, Trek VFI-10, Richard James FTS-50 and FTS-100 systems, Greishaber VIS, and others. Pack includes 10cc syringe with pressure-tested tubing and adapters, stainless steel cannula, and Teflon cannula. 0.2 micron filter included to maintain sterile field.

10 Packs per box, Sterile, Single-Use

Reorder Number TR8170-11

Irrigating Contact Lens

Ideal for use during silicone oil procedures, irrigating contact lens enables visualization of surgical field while ensuring proper irrigation of cornea and orbit. Available in flat or minus styles, single-use lens is convenient requiring no oil clean up.

Packaged individually, Sterile, Single Use

Reorder Number TR9840-11 (Flat Lens) / TR9841-11 (Minus Lens)



Viscous Fluid Cannula



Viscous Fluid Pack



Irrigating Contact Lens

