



Sonomed Escalon™

Intraocular Gases

SF6, C3F8, Gas Mixing Kits, and Accessories

SF6 Intraocular Gas [TR9088]

20 gram multi-dose cylinder, > 99.99% purity, 30 patient bracelets and cards.

C3F8 Intraocular Gas [TR9089]

20 gram multi-dose cylinder, >99.8% purity, 30 patient bracelets and cards.

Universal Gas Kits [TR9060-61]

Convenient kit for the controlled and accurate mixing and dispensing of gases. Includes 3-way stopcock for drawing of gas or air, 0.2 micron filters provide for sterilization of gas and/or air, and 60cc syringe pre-marked with mixture percentages to eliminate guesswork during prepping. Provided sterile, 10 per box.

Tornambe Eye Level [TR9049-11]

Convenient level aids in communicating and providing feedback to the patient to ensure proper positioning of the head following insertion of gas bubble.

Regulator [TR9085]

Cylinder Stand [TR9066]



-
- 01 C3F8 Gas Cylinder with Regulator
 - 02 Universal Gas Kit
 - 03 Tornambe Eye Level
-